


ResMed Ultra Mirage Nasal Interface Sizing Gauge

16984/1 03 04

STANDARD **LARGE** **1 inch** **SHALLOW** **SHALLOW-WIDE**

ALIGN WITH BOTTOM OF NOSE ALIGN WITH BOTTOM OF NOSE ALIGN WITH BOTTOM OF NOSE ALIGN WITH BOTTOM OF NOSE

Follow these two steps to choose the correct size:



1 Check the nostril width for clearance.

2 Check the nasal bridge width for clearance.

Ultra
MIRAGE™ Nasal Mask

FITTING TEMPLATE

RESMED

Customer Service 1-800-424-0737 (US only)
www.resmed.com

© Ultra Mirage and Mirage are trademarks of ResMed Ltd. © 2003 ResMed Ltd.

INSTRUCTIONS:

1. Print out
2. Measure inch marker to verify scale
3. Cut out grey areas
4. Choose size based on instruction printed on gauge