

## DeVilbiss IntelliPAP<sup>™</sup> AutoAdjust<sup>®</sup> CPAP System

# So **intelligent,** it automatically adjusts.

The new **DeVilbiss IntelliPAP AutoAdjust CPAP System** incorporates all of the design and performance advantages of the IntelliPAP Standard CPAP System in addition to a **wide range of customizable features**.

This system has been **designed with patients and providers in mind,** incorporating many of the recommendations solicited through research to optimize patient adherence. IntelliPAP's combination of **auto-adjusting pressure, comfort, education and adherence tracking with SmartCode™** helps to ensure patient compliance – all in a **highly efficient platform** that makes great business sense to you.

From its **reliable electronics and water resistant housing** to its durable Lexan<sup>®</sup> humidifier chamber and stainless steel heat transfer plate, IntelliPAP is **built for years of dependable performance and patient satisfaction.** Just what you would expect from DeVilbiss – the widely recognized leader in respiratory technology.

DeVilbiss IntelliPAP AutoAdjust CPAP Unit Model DV54D-HH

## Performance Features

- Exclusive DeVilbiss AutoAdjust technology
- Quietest CPAP Unit on the market ONLY 26 dBA!\*
- Integrated Heated Humidification
- Automatic Altitude Adjustment up to 9,000 feet
- Exclusive humidifier chamber design maximizes humidification and minimizes leaks
- Exceptionally small footprint, smallest integrated unit among top-selling CPAP
- Delay pressure ramp

Try our Serenity<sup>®</sup>, FlexSet<sup>®</sup> or FlexAire<sup>®</sup> masks – DeVilbiss offers a lightweight, comfortable and latex-free line of masks. Visit www.DeVilbissHealthcare.com for more information.

### Design Features

- Whisper-quiet operation 26 dBA\*
- Exceptionally small footprint fits easily into the home environment and carry-on luggage
- Discrete **carry bag** included that can be easily expanded to hold additional accessories
- Integrated heated humidification system easily docks into the bottom of the unit without requiring a separate hose or power cord
- Humidifier reservoir is easy to remove and replace, and is dishwasher safe
- Water ingress protection prevents water from humidifier damaging the CPAP – great for travelers
- Intuitive, **3-large buttons** control nearly all functions for simplified patient operation
- Backlit LCD and easy-to-follow patient menu
- Delay pressure ramp eases the machine to the prescribed pressure allowing the user to fall asleep at a lower pressure (0-45 minutes in 5 minute increments and .5 increments for pressure)
- Patient air tubing extends from unit back to reduce clutter and **enhance** visual appeal
- Accommodates various tubing lengths (6 and 10 feet)
- Eight-foot power cord length provides **greater positioning flexibility** - no bulky power brick on cord
- Clean, sophisticated design
- Large, slip-resistant pads on bottom of unit to minimize movement on night table during the night
- Auto on/off automatically starts the flow of air after breathing into the mask and stops the flow of air when mask is removed
- SmartCode<sup>™</sup> compliance tracking encodes usage data that can be communicated to a clinician and then translated
- Extra air inlet and fine particle filter included
- Industry leading, **3-year warranty**

#### **Technical Specifications**

Product Dimensions	6.4" H x 6.5" W x 8.4" D
Product Weight	4.45 lbs w/humidifier and PAP
Pressure Range	3-20 cm H20
Universal Power Supply Yes	
AC Cord Length	8 feet
Automatic Altitude Adjustment	Yes
DC Operable	Yes
ABOB Analog Output for PSG	Yes (pressure, patient flow, leak flow, patient event markers, est. tidal volume)

Automatic Leak Compensation	Yes
Auto On/Off	Yes (can be enabled or disabled)
Motor Hour Meter	Yes
Therapy Compliance Meter	Yes
Warranty	3-year warranty on IntelliPAP 2-year warranty on humidification system
DeVilbiss Smart Link Compatible	Yes
SmartCode Feature	Yes



DeVilbiss Healthcare100 DeVilbiss DriveSomerset, PA 15501800.338.1988814.443.4881www.DeVilbissHealthcare.com

Patient Education and Marketing Support

- User video includes an introduction to CPAP therapy and tips for optimizing therapy and how to care for and use the IntelliPAP system
- Fully-illustrated patient guide demonstrates how easy it is to use IntelliPAP, while providing a comprehensive overview of sleep apnea and its treatment

#### AutoAdjust Features

The IntelliPAP AutoAdjust unit can be operated in **both conventional CPAP and AutoAdjust modes.** AutoAdjust features include:

- Automatic pressure adjustment based on event density and type
- Adjustable apnea/hypopnea detection
- Adjustable upper and lower pressure limits
- Adjustable delay setting for up to four hours ideal for comparing patient response to AutoAdjust and standard CPAP therapy
- Automatic detection of mask conditions and ultra-sensitive snoring detection



DeVilbiss IntelliPAP AutoAdjust CPAP Unit

\* Tested per ISO 17510

©2007 DeVilbiss Healthcare 11.07 DeVilbiss Healthcare is a division of Sunrise Medical HHG, Inc. IntelliPAP, Serenity, FlexSet and FlexAire are registered trademarks of DeVilbiss Healthcare. LEXAN is a registered trademark of GE Plastics. LT-2009 Rev. A